APPLICATION NOTE

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Table 1. Preparation of Calibration Standards

Calibration Standard (% Phosphorylated)	μL Peptide Stock	μL Phosphopeptide Stock
0	100	0
12.5	87.5	12.5
25	75	25
50	50	50
100	0	100

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

Table 2. IMAP FP Settings in SoftMax Pro for SpectraMax M5 Microplate Reader

Read type	Endpoint	
Read mode	Fluorescence polarization	
	FAM	TAMRA
Wavelengths	Ex 485 nm Em 525 nm Cutoff 515 nm	Ex 530 nm Em 590 nm Cutoff 570 nm
Sensitivity	Readings: 100 PMT: high or medium	
Automix	Off	
AutoCalibrate	On	
Assay plate type	384-well Costar black	
Wells to read	[Determined by user]	
Settling time	Delay: 100 ms	
AutoRead	Off	

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IMAP-FP kinase assay on SpectraMax M5 (figure 2)

Lck kinase dilution curve with FAM-p34cdc2-derived peptide from 0.04 to 3.5 units/mL enzyme, read on SpectraMax M5 (4-parameter curve fit). Error bars here and on all subsequent graphs are standard deviations.

IMAP-FP kinase assay on Analyst HT (figure 3)

Lck kinase dilution curve with FAM-p34cdc2-derived peptide from 0.04 to 3.5 units/mL enzyme, read on Analyst HT (4-parameter curve fit).

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


IMAP with FAM-labeled peptide (figure 4)


Calibration standard curve with FAM-p34cdc2-derived peptide and phosphopeptide, read on SpectraMax M5 (quadratic curve fit).

Lck kinase assay with TAMRA-labeled peptide on SpectraMax M5 (figure 5)

Lck kinase dilution curve with TAMRA-p34cdc2-derived peptide, read on SpectraMax M5 (4-parameter curve fit).

Lck kinase assay with TAMRA-labeled peptide on Analyst HT (figure 6)

Lck kinase dilution curve with TAMRA-p34cdc2-derived peptide, read on Analyst HT (4-parameter curve fit).

IMAP with TAMRA-labeled peptide (figure 7)

Calibration standard curve with TAMRA-p34cdc2-derived peptide and phosphopeptide, read on SpectraMax M5 (linear curve fit).

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